

|  |  |          |
|--|--|----------|
| Harris County Institute of Forensic Sciences |  |          |
| Section: Toxicology                          | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  | Rev.: 21 |

1.0 Purpose

1.1 This document lists the major analytes anticipated to produce positive responses in the routine screens applied in the Forensic Toxicology Laboratory.

2.0 Scope

2.1 This procedure is applicable to the Forensic Toxicology Laboratory

2.2 The routine screening analytes are evaluated annually. For example:

2.2.1 Analytes may be added if emerging substances are detected by other techniques or in other laboratories.

2.2.2 Substances may be removed if they are rarely detected in casework and/or determined to be of minor significance.

3.0 Definitions and Abbreviations

3.1 No method-specific or non-standard terms are used in this procedure.

4.0 Materials

4.1 Not Applicable

5.0 Procedure

5.1 Tables

5.1.1 Volatiles

| Analyte     | Established Cut-off |
|-------------|---------------------|
| Ethanol     | 0.010 g/100 mL      |
| Acetone     | 0.010 g/100 mL      |
| Methanol    | 0.010 g/100 mL      |
| Isopropanol | 0.010 g/100 mL      |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

5.1.2 ELISA (cross-reactivity with structurally similar compounds in **blood** from ELISA inserts, Immunalysis Corporation)

| Analyte  | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 10 ng/mL THCA | Cross-Reactivity (%) |
|--|---|------------------------------------|---|----------------------|
| <b>(-)-11-nor-9-Carboxy-<math>\Delta</math>9-THC</b> | 10  | N/A                                | 10  | 100                  |
| 11-hydroxy- $\Delta$ 9-THC                           |   |                                    | 200                                       | 5                    |
| (-) $\Delta$ 8-THC                                   |   |                                    | 300                                       | 3.33                 |
| (-) $\Delta$ 9-THC                                   |   |                                    | 400                                       | 2.5                  |
| Cannabinol   |   |                                    | 1,000                                     | 1                    |
| Cannabidiol  |   |                                    | 300,000                                   | 0.003                |
| (+) 11-carboxylic acid gluc                          |   |                                    | 1,000                                     | 1                    |

| Analyte                         | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 10 ng/mL Morphine | Cross-Reactivity (%) |
|---------------------------------|---|------------------------------------|---|----------------------|
| <b>Morphine</b>                 | 20  | 20                                 | 10  | 100.00%              |
| 6-Acetyl Codeine                |   |                                    | 15  | 66.70%               |
| 6-Acetyl Morphine               |   |                                    | 40  | 25.00%               |
| Codeine                         |   |                                    | 5   | 200.00%              |
| Desomorphine                    |   |                                    | 18  | 55.60%               |
| Dihydrocodeine                  |   |                                    | 20  | 50.00%               |
| Ethylmorphine                   |   |                                    | 5   | 200.00%              |
| Heroin                          |   |                                    | 100   | 10.00%               |
| Hydrocodone                     |   |                                    | 10  | 100.00%              |
| Hydromorphone                   |   |                                    | 20  | 50.00%               |
| Levorphanol                     |   |                                    | 65  | 15.40%               |
| Morphine-3 $\beta$ -Glucuronide |   |                                    | 12.5  | 80.00%               |
| Morphine-6 $\beta$ -Glucuronide |   |                                    | 75,000  | 0.01%                |
| Nalorphine                      |   |                                    | 2,200   | 0.50%                |
| Naloxone                        |   |                                    | 100,000                                       | 0.01%                |
| Norcodeine                      |   |                                    | 1,400   | 0.70%                |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                |  |  |         |        |
|----------------|--|--|---------|--------|
| Normorphine    |  |  | 15,000  | 0.10%  |
| Noroxycodone   |  |  | 80,000  | 0.01%  |
| Noroxymorphone |  |  | 100,000 | <0.01% |
| Oxycodone      |  |  | 1,500   | 0.70%  |
| Oxymorphone    |  |  | 3,000   | 0.30%  |

| Analyte                                | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 5 ng/mL PCP | Cross-Reactivity (%) |
|--|---|------------------------------------|---|----------------------|
| <b>Phencyclidine (PCP)</b>             | 10  | 20                                 | 5                                       | 100                  |
| 1-[1-(2-thienyl)cyclohexyl]-piperidine |   |                                    | 10                                      | 50                   |
| 1-[1-(2-thienyl)cyclohexyl]-morpholine |   |                                    | 50                                      | 10                   |
| 1-(1-phenylcyclohexyl)pyrrolidine      |   |                                    | 50                                      | 10                   |
| 1-(1-phenylcyclohexyl)morpholine       |   |                                    | 60                                      | 8                    |

| Analyte                | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 20 ng/mL S(+)-Amphetamine | Cross-Reactivity (%) |
|------------------------|---|------------------------------------|---|----------------------|
| <b>(+) Amphetamine</b> | 20  | 40                                 | 20  | 100                  |
| (-)Amphetamine         |   |                                    | 1500  | 1.3                  |
| (±)Amphetamine         |   |                                    | 50  | 40                   |
| Diphenhydramine        |   |                                    | 100,000   | 0.02                 |
| (+)Ephedrine           |   |                                    | 100,000   | 0.02                 |
| (-)Ephedrine           |   |                                    | 100,000   | 0.02                 |
| Fenfluramine           |   |                                    | 100,000   | 0.02                 |
| 4-Fluoroamphetamine    |   |                                    | 100   | 20                   |
| MBDB                   |   |                                    | 20,000  | 0.1                  |
| MDA                    |   |                                    | 15  | 133.3                |
| MDEA                   |   |                                    | 3,500   | 0.6                  |
| MDMA                   |   |                                    | 5,000   | 0.4                  |
| (+)Methamphetamine     |   |                                    | 20,000  | 0.1                  |
| (-)Methamphetamine     |   |                                    | 100,000   | 0.02                 |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                     |  |  |         |      |
|---------------------|--|--|---------|------|
| Norpseudoephedrine  |  |  | 20,000  | 0.1  |
| Phentermine         |  |  | 150     | 13.3 |
| Phenylephrine       |  |  | 100,000 | 0.02 |
| Phenylpropanolamine |  |  | 70,000  | 0.03 |
| PMA                 |  |  | 10      | 200  |
| (+)Pseudoephedrine  |  |  | 100,000 | 0.02 |
| (-)Pseudoephedrine  |  |  | 100,000 | 0.02 |
| Tyramine            |  |  | 6,000   | 0.3  |

| Analyte                   | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 50 ng/mL Benzoyllecgonine | Cross-Reactivity (%) |
|---------------------------|---|------------------------------------|---|----------------------|
| <b>Benzoyllecgonine</b>   | 20  | 20                                 | 50  | 100                  |
| Cocaethylene              |   |                                    | 50  | 100                  |
| Cocaine                   |   |                                    | 60  | 83                   |
| Ecgonine                  |   |                                    | 100,000   | 0.1                  |
| Ecgonine Ethyl Ester      |   |                                    | 55,000  | 0.1                  |
| Ecgonine Methyl Ester     |   |                                    | 25,000  | 0.2                  |
| m-Hydroxybenzoyllecgonine |   |                                    | 100   | 50                   |
| Norcocaine                |   |                                    | 55,000  | 0.1                  |

| Analyte          | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 10 ng/mL Oxycodone | Cross-Reactivity (%) |
|------------------|---|------------------------------------|--|----------------------|
| <b>Oxycodone</b> | 5   | 5                                  | 10   | 100                  |
| 6-Acetylmorphine |   |                                    | 8,000  | 0.13                 |
| Buprenorphine    |   |                                    | 100,000  | <0.01                |
| Codeine          |   |                                    | 8,000  | 0.13                 |
| Dihydrocodeine   |   |                                    | 7,000  | 0.14                 |
| Heroin           |   |                                    | 10,000   | 0.1                  |
| Hydrocodone      |   |                                    | 400  | 2.5                  |
| Hydromorphone    |   |                                    | 500  | 2.0                  |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                                    |  |  |         |       |
|------------------------------------|--|--|---------|-------|
| Levorphanol                        |  |  | 5,000   | 0.2   |
| Meperidine                         |  |  | 100,000 | <0.01 |
| Morphine                           |  |  | 10,000  | 0.1   |
| Morphine-3-Glucuronide             |  |  | 10,000  | 0.1   |
| Nalbuphine                         |  |  | 100,000 | <0.01 |
| Nalorphine                         |  |  | 100,000 | <0.01 |
| Naloxone                           |  |  | 10,000  | 0.1   |
| Naloxone-3-Glucuronide             |  |  | 100,000 | <0.01 |
| Naltrexone                         |  |  | 100,000 | <0.01 |
| Norcodeine                         |  |  | 100,000 | <0.01 |
| Normorphine                        |  |  | 100,000 | <0.01 |
| Noroxycodone                       |  |  | 200     | 5.0   |
| Noroxymorphone                     |  |  | 500     | 2.0   |
| Oxymorphone                        |  |  | 10      | 100   |
| Oxymorphone-3 $\beta$ -Glucuronide |  |  | 400     | 2.5   |
| Tramadol                           |  |  | 100,000 | N.D.  |

| Analyte                    | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 20 ng/mL S (+)-Methamphetamine | Cross-Reactivity (%) |
|----------------------------|---|------------------------------------|--|----------------------|
| <b>(+) Methamphetamine</b> | 20  | 20                                 | 20   | 100                  |
| (-)Amphetamine             |   |                                    | 100,000  | 0.02                 |
| (+)Amphetamine             |   |                                    | 25,000   | 0.1                  |
| Diphenhydramine            |   |                                    | 100,000  | 0.02                 |
| (-)Ephedrine               |   |                                    | 1,000  | 2                    |
| (+)Ephedrine               |   |                                    | 100,000  | 0.02                 |
| Fenfluramine               |   |                                    | 5,000  | 0.4                  |
| MBDB                       |   |                                    | 2,000  | 1                    |
| MDA                        |   |                                    | 10,000   | 0.2                  |
| MDMA                       |   |                                    | 20   | 100                  |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                     |  |  |         |      |
|---------------------|--|--|---------|------|
| MDEA                |  |  | 40      | 50   |
| (-) Methamphetamine |  |  | 3,000   | 0.7  |
| Norpseudoephedrine  |  |  | 100,000 | 0.02 |
| Phentermine         |  |  | 100,000 | 0.02 |
| Phenylephrine       |  |  | 100,000 | 0.02 |
| Phenylpropanolamine |  |  | 100,000 | 0.02 |
| PMA                 |  |  | 5,000   | 0.4  |
| (+)Pseudoephedrine  |  |  | 50,000  | 0.04 |
| (-)Pseudoephedrine  |  |  | 100,000 | 0.02 |
| Tyramine            |  |  | 100,000 | 0.02 |

| Analyte                  | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 50 ng/mL Oxazepam | Cross-Reactivity (%) |
|--------------------------|---|------------------------------------|---|----------------------|
| <b>Oxazepam</b>          | 10  | 20                                 | 50  | 100.0                |
| Alprazolam               |   |                                    | 18  | 277.8                |
| 7-Aminoclonazepam        |   |                                    | 180   | 27.8                 |
| Bromazepam               |   |                                    | 55  | 90.9                 |
| Chlordiazepoxide         |   |                                    | 275   | 18.2                 |
| Clobazam                 |   |                                    | 45  | 111.1                |
| Clonazepam               |   |                                    | 65  | 76.9                 |
| Clorazepate              |   |                                    | 30  | 166.7                |
| Desalkylflurazepam       |   |                                    | 30  | 166.7                |
| N-Desmethylflunitrazepam |   |                                    | 40  | 125.0                |
| Diazepam                 |   |                                    | 25  | 200.0                |
| Estazolam                |   |                                    | 25  | 200.0                |
| Etizolam                 |   |                                    | 40  | 125.0                |
| Flubromazepam            |   |                                    | 25  | 200.0                |
| Flubromazolam            |   |                                    | 20  | 250.0                |
| Flurazepam               |   |                                    | 30  | 166.7                |
| Flunitrazepam            |   |                                    | 40  | 125.0                |
| Flumazenil               |   |                                    | 100,000                                       | N.D.                 |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                              |  |  |        |         |
|------------------------------|--|--|--------|---------|
| $\alpha$ -Hydroxy-Alprazolam |  |  | 25     | 200.0   |
| 2-Hydroxyethylflurazepam     |  |  | 30     | 16670.0 |
| Hydroxy-Triazolam            |  |  | 40     | 125.0   |
| Lorazepam Glucuronide        |  |  | 50,000 | 0.1     |
| Lorazepam                    |  |  | 80     | 62.5    |
| Lormetazepam                 |  |  | 50     | 100.0   |
| Oxazepam Glucuronide         |  |  | 10,000 | 0.5     |
| Medazepam                    |  |  | 85     | 58.8    |
| Midazolam                    |  |  | 45     | 111.1   |
| Nimetazepam                  |  |  | 30     | 166.7   |
| Nitrazepam                   |  |  | 50     | 100.0   |
| Nordiazepam                  |  |  | 30     | 166.7   |
| Prazepam                     |  |  | 35     | 142.9   |
| Temazepam                    |  |  | 35     | 142.9   |
| Triazolam                    |  |  | 35     | 142.9   |

| Analyte       | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 50 ng/mL Methadone | Cross-Reactivity (%) |
|---------------|---|------------------------------------|--|----------------------|
| (±) Methadone | 20  | 20                                 | 50   | 100                  |
| Methadol      |   |                                    | 600  | 8.3                  |
| EDDP          |   |                                    | 50,000   | 0.1                  |
| EMDP          |   |                                    | 100,000  | <0.1                 |
| LAAM          |   |                                    | 4,000  | 1.3                  |

| Analyte            | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 1 ng Fentanyl | Cross-Reactivity (%) |
|--------------------|---|------------------------------------|---|----------------------|
| <b>Fentanyl</b>    | 1   | 2                                  | 1.0                                       | 100.0                |
| Acetyl fentanyl    |   |                                    | 1.8                                       | 55.6                 |
| Acetyl Norfentanyl |   |                                    | 1,000,000                                 | <0.0001              |
| Acrylfentanyl      |   |                                    | 2.0                                       | 50.0                 |

|  |  |  |          |
|--|--|--|----------|
| Harris County Institute of Forensic Sciences |  |  |          |
| Section: Toxicology                          |  | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         |  | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  |  | Rev.: 21 |

|                                    |  |  |           |         |
|------------------------------------|--|--|-----------|---------|
| Alfentanil                         |  |  | 1,000,000 | <0.0001 |
| Beta-Hydroxy fentanyl              |  |  | 100       | 1.0     |
| Butyryl fentanyl                   |  |  | 1.0       | 100     |
| Carfentanil                        |  |  | 1000      | 0.1     |
| Cis-d,l-3-methylfentanyl           |  |  | 50        | 2.0     |
| Despropionylfentanyl (4-ANPP)      |  |  | 50        | 2.0     |
| FIBF                               |  |  | 3.0       | 33.3    |
| Furanyl fentanyl                   |  |  | 2.0       | 50.0    |
| Isobutyryl fentanyl                |  |  | 2.0       | 50.0    |
| meta-Fluorobutyryl fentanyl        |  |  | 2.0       | 50.0    |
| meta-Fluorofentanyl                |  |  | 2.0       | 50.0    |
| Methoxyacetyl fentanyl             |  |  | 1.5       | 66.7    |
| 4-methoxy-butyryl fentanyl         |  |  | 1.8       | 55.6    |
| meta-Chlorophenylpiperazine (mCPP) |  |  | 1,000,000 | 0.0001  |
| Norfentanyl                        |  |  | 1,000,000 | 0.0001  |
| Ortho-Fluorobutyryl fentanyl       |  |  | 2.0       | 50.0    |
| Ortho-Fluorofentanyl               |  |  | 2.0       | 50.0    |
| para-Chloroisobutyryl fentanyl     |  |  | 4.0       | 25.0    |
| para-Fluorobutyryl fentanyl        |  |  | 1.75      | 57.1    |
| Remifentanil                       |  |  | 1,000,000 | <0.0001 |
| Sufentanil                         |  |  | 1000      | 0.1     |
| Trans-d,l-3-methylfentanyl         |  |  | 5.0       | 20.0    |
| U-47700                            |  |  | 25,000    | 0.004   |
| Valeryl Fentanyl                   |  |  | 1.5       | 50.0    |
| W-18                               |  |  | 1,000,000 | <0.0001 |

| Analyte             | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | Approx. ng/mL equivalent to 1000 ng Carisoprodol | Cross-Reactivity (%) |
|---------------------|---|------------------------------------|--|----------------------|
| <b>Carisoprodol</b> | 500                                       | 500                                | 1000   | 100                  |
| Meprobamate         |   |                                    | 850  | 118                  |



|  |  |  |
|--|--|--|
| Harris County Institute of Forensic Sciences |  |  |
| Section: Toxicology                          | Approved by: Director, Forensic Toxicology |  |
| Document Type: QA/QC                         | Procedure No.: TOX07.1027                  |  |
| Title: Toxicology Routine Screening Analytes | Rev.: 21                                   |  |

| Analyte       | Established Cut-off, ng/mL (Blood, Urine) | Established Cut-off, ng/mL (Liver) | No Cross-Reactivity Data |
|---------------|---|------------------------------------|--------------------------|
| Acetaminophen | 1000                                      | N/A                                |                          |

### 5.1.3 GC/MS Screens

| Analyte                            | Concentration in Control | Administrative Cut-off |
|------------------------------------|--------------------------|------------------------|
| <b>Standard Acidic Drug Screen</b> |                          |                        |
| Butalbital                         | 10 mg/L                  | 2 mg/L                 |
| Carbamazepine                      | 2 mg/L                   | 1 mg/L                 |
| Carisoprodol                       | 10 mg/L                  | 2 mg/L                 |
| Ibuprofen                          | 25 mg/L                  | 2 mg/L                 |
| Levetiracetam                      | 10 mg/L                  | 2 mg/L                 |
| Meprobamate                        | 25 mg/L                  | 2 mg/L                 |
| Naproxen                           | 25 mg/L                  | 2 mg/L                 |
| Pentobarbital                      | 10 mg/L                  | 2 mg/L                 |
| Phenobarbital                      | 10 mg/L                  | 2 mg/L                 |
| Phenytoin                          | 10 mg/L                  | 2 mg/L                 |
| Secobarbital                       | 10 mg/L                  | 2 mg/L                 |
| Topiramate                         | 10 mg/L                  | 2 mg/L                 |
| Valproic Acid                      | 10 mg/L                  | 2 mg/L                 |
| <b>Standard Basic Drug Screen</b>  |                          |                        |
| Alprazolam                         | 0.25 mg/L                | N/A                    |
| Amitriptyline                      | 0.25 mg/L                | 0.03 mg/L              |
| Carisoprodol                       | 10 mg/L                  | 2 mg/L                 |
| Chlorpheniramine                   | 0.10 mg/L                | 0.03 mg/L              |
| Citalopram                         | 0.25 mg/L                | 0.03 mg/L              |
| Cocaine                            | 0.25 mg/L                | 0.03 mg/L              |
| Codeine                            | 0.25 mg/L                | 0.03 mg/L              |
| Cyclobenzaprine                    | 0.25 mg/L                | 0.03 mg/L              |
| Dextromethorphan (Methorphan)      | 0.10 mg/L                | 0.03 mg/L              |
| Diltiazem                          | 0.25 mg/L                | 0.03 mg/L              |
| Diazepam                           | 0.25 mg/L                | 0.03 mg/L              |

|  |  |          |
|--|--|----------|
| Harris County Institute of Forensic Sciences |  |          |
| Section: Toxicology                          | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  | Rev.: 21 |

|                           |           |           |
|---------------------------|-----------|-----------|
| Diphenhydramine           | 0.25 mg/L | 0.03 mg/L |
| Doxylamine                | 0.25 mg/L | 0.03 mg/L |
| Fentanyl                  | 0.02 mg/L | N/A       |
| Hydrocodone               | 0.10 mg/L | 0.03 mg/L |
| Levamisole (Tetramisole)  | 0.25 mg/L | 0.03 mg/L |
| Lidocaine                 | 0.25 mg/L | 0.03 mg/L |
| Meprobamate               | 25 mg/L   | 2 mg/L    |
| Methadone                 | 0.25 mg/L | 0.03 mg/L |
| Methamphetamine           | 0.25 mg/L | 0.03 mg/L |
| Mirtazapine               | 0.25 mg/L | 0.03 mg/L |
| Nordiazepam               | 0.25 mg/L | 0.03 mg/L |
| Norsertaline              | 0.25 mg/L | 0.03 mg/L |
| Nortriptyline             | 0.25 mg/L | 0.03 mg/L |
| Oxycodone                 | 0.10 mg/L | N/A       |
| Phencyclidine             | 0.25 mg/L | 0.03 mg/L |
| Promethazine              | 0.25 mg/L | 0.03 mg/L |
| Propoxyphene              | 0.25 mg/L | 0.03 mg/L |
| Sertraline                | 0.25 mg/L | 0.03 mg/L |
| Tramadol                  | 0.25 mg/L | 0.03 mg/L |
| Trazodone                 | 0.25 mg/L | 0.03 mg/L |
| Venlafaxine               | 0.25 mg/L | 0.03 mg/L |
| Verapamil                 | 0.25 mg/L | 0.03 mg/L |
| Zolpidem                  | 0.25 mg/L | 0.03 mg/L |
| <b>Other Drug Screens</b> |           |           |
| GHB                       | 5 mg/L    | 5 mg/L    |

5.1.4 TOF Screen:

| Compound             | Concentration |                                    |
|----------------------|---------------|------------------------------------|
| 25I-NBOMe            | 10 ng/mL      |                                    |
| 7-Aminoflunitrazepam | 5 ng/mL       |                                    |
| Brompheniramine      | 10 ng/mL      |                                    |
| Buprenorphine        | 1 ng/mL       | Blood only; not approved for urine |
| Butylone             | 10 ng/mL      |                                    |
| Carbamazepine        | 1000 ng/mL    |                                    |

|  |  |  |
|--|--|--|
| Harris County Institute of Forensic Sciences |  |  |
| Section: Toxicology                          | Approved by: Director, Forensic Toxicology |  |
| Document Type: QA/QC                         | Procedure No.: TOX07.1027                  |  |
| Title: Toxicology Routine Screening Analytes | Rev.: 21                                   |  |

|                  |            |                                    |
|------------------|------------|------------------------------------|
| Chlordiazepoxide | 10 ng/mL   |                                    |
| Chlorpheniramine | 10 ng/mL   |                                    |
| Dextromethorphan | 10 ng/mL   |                                    |
| Diphenhydramine  | 10 ng/mL   |                                    |
| Ethylone         | 10 ng/mL   |                                    |
| Ketamine         | 10 ng/mL   |                                    |
| Lorazepam        | 5 ng/mL    |                                    |
| MDPV             | 10 ng/mL   |                                    |
| Mephedrone       | 10 ng/mL   |                                    |
| Methedrone       | 50 ng/mL   |                                    |
| Methylone        | 10 ng/mL   |                                    |
| Methylphenidate  | 50 ng/mL   |                                    |
| Midazolam        | 10 ng/mL   |                                    |
| MPPP             | 10 ng/mL   |                                    |
| Naphyrone        | 10 ng/mL   |                                    |
| Oxcarbazepine    | 1000 ng/mL |                                    |
| Oxycodone        | 10 ng/mL   |                                    |
| Oxymorphone      | 10 ng/mL   |                                    |
| Promethazine     | 10 ng/mL   |                                    |
| Quetiapine       | 10 ng/mL   |                                    |
| Ritalinic Acid   | 100 ng/mL  |                                    |
| TFMPP            | 50 ng/mL   |                                    |
| TMA-2            | 50 ng/mL   |                                    |
| Topiramate       | 50 ng/mL   |                                    |
| Tramadol         | 10 ng/mL   |                                    |
| Zaleplon         | 10 ng/mL   |                                    |
| Zolpidem         | 10 ng/mL   | Blood only; not approved for urine |
| Zopiclone        | 10 ng/mL   |                                    |

6.0 Data Analysis and Documentation

6.1 Not Applicable

7.0 Acceptance Criteria

7.1 Not Applicable

|  |  |          |
|--|--|----------|
| Harris County Institute of Forensic Sciences |  |          |
| Section: Toxicology                          | Approved by: Director, Forensic Toxicology |          |
| Document Type: QA/QC                         | Procedure No.: TOX07.1027                  |          |
| Title: Toxicology Routine Screening Analytes |  | Rev.: 21 |

8.0 References

8.1 ELISA inserts, Immualysis Corporation

9.0 Attachment or Appendices

9.1 Not applicable